



QUALITY CLINICAL CHEMISTRY RANGE FOR BRIDGING THE HEALTHIER FUTURE

ALLCARE Clinical Chemistry Reagents are designed to enhance user experience by facilitating ease of use and quality results at affordable rates.

QUICK FEATURES:

- ◆ Ready to use Liquid stable Reagents
- ◆ Accurate results and high linearity
- ◆ Available in small, medium and bulk pack
- ◆ Reagents in Monoval for sensitive parameters
- ◆ Kits provided with standard / calibrator
- ◆ Easy program and user friendly procedure





CLINICAL CHEMISTRY PRODUCT LIST

Product Code	Product Name	Pack Size	Method	Type
ABI01001	ADA	1 x 15 ml	ENZYMATIC KINETIC	LIQ. STABLE 2 REAGENT
ABI01002		2 x 15ml		
ABI02001	ALBUMIN	1 x 100 ml	BCG	LIQ. STABLE 1 REAGENT
ABI02002		2 x 100 ml		
ABI03001	ALKALINE PHOSPHATASE (SR)	25 x 1.1 ml	PNPP/AMP KINETIC	LIQ. STABLE 1 REAGENT
ABI03003		50 x 1.1 ml		
ABI03002		4 x 25 ml		
ABI40001	ALKALINE PHOSPHATASE (DR)	2 x 24 ml	PNPP/AMP KINETIC	LIQ. STABLE 2 REAGENT
ABI40002		4 x 25 ml		
ABI04001	AMYLASE	1 x 10 ml	GAL-G2-CNP	LIQ. STABLE 1 REAGENT
ABI04002		2 x 10 ml		
ABI05001	BILIRUBIN (T & D) SR	4 x 25 ml	BERGH & MULLER	LIQ. STABLE 1 REAGENT
ABI05002		4 x 50 ml		
ABI39001	BILIRUBIN (T & D)	4 x 50 ml	JENDRASSIK-GROF	LIQ. STABLE 2 REAGENT
ABI39002	BILIRUBIN TOTAL	4 x 50 ml		
ABI39003	BILIRUBIN DIRECT	4 x 50 ml		
ABI06001	CALCIUM	25 x 1.1 ml	ARSENazo III	LIQ. STABLE 1 REAGENT
ABI06003		50 x 1.1 ml		
ABI06002		4 x 25 ml		
ABI07001	TOTAL CHLOLESTEROL	4 x 25 ml	CHOD/POD	LIQ. STABLE 1 REAGENT
ABI07002		4 x 50 ml		
ABI08001	CHLORIDE	25 x 1.1 ml	MTC	LIQ. STABLE 1 REAGENT
ABI08003		50 x 1.1 ml		
ABI08002		4 x 50 ml		
ABI09001	CREATININE (SR)	2 x 50 ml	ALK-PICRATE / JEFFE	LIQ. STABLE 1 REAGENT
ABI09002		4 x 50 ml		
ABI38001	CREATININE (DR)	2 x 50 ml	ALK-PICRATE / JEFFE	LIQ. STABLE 2 REAGENT
ABI38002		4 x 50 ml		
ABI10001	CRP	1 x 50 ml	LATEX TURBIDIMETRY	LIQ. STABLE 2 REAGENT
ABI10002		2 x 50ml		
ABI11001	CK-NAC	1 x 10 ml	NAC-ACTIVATED	LIQ. STABLE 2 REAGENT
ABI11002		2 x 10 ml		
ABI12001	CK-MB	1 x 10 ml	IMMUNO INHIBITION	LIQ. STABLE 2 REAGENT
ABI12002		2 x 10 ml		
ABI13001	D-Dimmer with calibrator	1 x 40 ml	LATEX TURBIDIMETRY	LIQ. STABLE 2 REAGENT
ABI14001	GLUCOSE	6 x 100 ml	GOD/POD	LIQ. STABLE 1 REAGENT
ABI14002		2 x 500 ml		
ABI15001	HDL CHOLESTROL	2 x 20 ml	SDT	LIQ. STABLE 2 REAGENT
ABI15002		2 x 40 ml		

CLINICAL CHEMISTRY PRODUCT LIST

Product Code	Product Name	Pack Size	Method	Type
ABI16001	LDH	1 x 20 ml	IFFC KINETIC	LIQ. STABLE 2 REAGENT
ABI16002		2 x 30 ml		
ABI17001	LDL CHOLESTROL	2 x 20 ml	SDT	LIQ. STABLE 2 REAGENT
ABI17002		2 x 40 ml		
ABI18001	MAGNESIUM	25 x 1.1 ml	XYLIDYL BLUE	LIQ. STABLE 1 REAGENT
ABI18002		4 x 25 ml		
ABI19001	POTASSIUM (Mono)	25 x 1.1 ml	COLORIMETRIC	LIQ. STABLE 1 REAGENT
ABI19002		50 x 1.1 ml		
ABI20001	SGOT / AST	4 x 25 ml	IFCC KINETIC	LIQ. STABLE 2 REAGENT
ABI20002		4 x 50 ml		
ABI21001	SGPT / ALT	4 x 25 ml	IFCC KINETIC	LIQ. STABLE 2 REAGENT
ABI21002		4 x 50 ml		
ABI22001	TOTAL PROTEIN	2 x 100 ml	BIURET	LIQ. STABLE 1 REAGENT
ABI23001	TRIGLYCERIDES	4 x 25 ml	GPO/POD	LIQ. STABLE 1 REAGENT
ABI23002		4 x 50 ml		
ABI24001	UREA BIRTHELOT	4 x 50 ml	BERTHELOT	LIQ. STABLE 1 REAGENT
ABI25001	UREA UV	4 x 25 ml	UV KINETIC	LIQ. STABLE 2 REAGENT
ABI25002		4 x 50 ml		
ABI26001	URIC ACID	50 x 1.1 ml	URICASE/POD	LIQ. STABLE 1 REAGENT
ABI26002		25 x 1.1 ml		
ABI26003		4 x 25 ml		
ABI27001	SODIUM (Mono)	25 x 1.1 ml	COLORIMETRIC	LIQ. STABLE 1 REAGENT
ABI27002		50 x 1.1 ml		
ABI28001	LIPASE	1 x 10 ml	METHYL RESORUFIN	LIQ. STABLE 2 REAGENT
ABI28002		2 x 10 ml		
ABI29001	AMMONIA*	1 x 10 ml	UV KINETIC	LIQ. STABLE 2 REAGENT
ABI29002		2 x 10 ml		
ABI30001	PHOSPHORUS	25 x 1.1 ml	UV METHOD	LIQ. STABLE 1 REAGENT
ABI30002		50 x 1.1 ml		
ABI31001	G6PDH*	10 Tests		Powder Form
ABI32001	GAMMA GT	1 x 10 ml	IFCC KINETIC	LIQ. STABLE 2 REAGENT
ABI32002		2 x 10 ml		
ABI33001	CHOLINESTERASE	1 x 10 ml	BUTYRYLTHIOCOLINE	LIQ. STABLE 2 REAGENT
ABI33002		2 x 10 ml		
ABI34001	RF *	1 x 50 ml	LATEX TURBIDIMETRY	LIQ. STABLE 2 REAGENT
ABI35001	ASO*	1 x 50 ml	LATEX TURBIDIMETRY	LIQ. STABLE 2 REAGENT
ABI36001	HOMOCYSTEINE*			
ABI37001	HbA1c*			

*COMING SOON

QUALITY AND RELIABLE INSTRUMENT RANGE



ALLCARE AC200 Semi Auto Biochemistry Analyzer

- ◆ Reliable, Robust and Easy to Operate
- ◆ Colour Graphical Display
- ◆ 48 Direct access test keys on board
- ◆ 248 Tests can be programmed
- ◆ 7 Interference Filter range 340-670nm with 3 Optional Positions
- ◆ Test Mode: Endpoint/Two Point /Kinetic/ Absorbance

ALLCARE AIFA 100 IFA POCT Analyzer

- ◆ Compact, Quick and Easy to Operate
- ◆ Accurate & Fully Quantitative Results
- ◆ 7" True Colour LCD with Touch Screen
- ◆ Upto 100 Tests/Hr (Quick Mode)
- ◆ Stores upto 1,00,000 Tests results



ALLCARE AC9000 Series Electrolyte Analyzer

- ◆ Based on Ion Selective Electrode (ISE) Method
- ◆ Parameters: Na⁺, K⁺, Cl⁻, iCa⁺⁺ & pH
- ◆ Colour display with Touch Screen
- ◆ Long Life maintenance free electrode
- ◆ Alarm for Reagent used up, Air bubble, Low sample and Waste full

H8 Fully Automated HbA1c Analyzer

- ◆ Based on Gold Standard HPLC Method
- ◆ IFCC & NGSP approved
- ◆ Results within 130 Seconds
- ◆ 10 Position Auto sampler
- ◆ HbA1c results are excluding interference of variant and unstable factors
- ◆ Graph based results for HbA1a, HbA1b, HbF, La1c, HbA1c, HbA0, Avg Glucose, etc



Hemax 330 Hematology Analyzer

- ◆ 3 Part Diff with Double Chamber
- ◆ 20 Parameters + 3 Histograms (WBC, RBC & PLT)
- ◆ Throughput: 60 Tests per hour
- ◆ Store upto 50,000 test results
- ◆ LCD Display with Touch screen
- ◆ USB Port to attach Mouse and Keyboard

